

5214 RESCO MULTILAB LABORATORY BENCHTOP SYSTEM

This section deals with the supply and installation of Resco proprietary laboratory benchtop system.

1. GENERAL

1.1 MANUFACTURER'S DOCUMENTS

Copies of current literature are available by contacting **RESCO** Ltd.

1.2 SAMPLES

Submit on request samples of panels and components for approval of colour and finishes.

1.3 QUALIFICATIONS

Workers to be experienced, competent and familiar with the materials and techniques specified. Installation to be carried out by workers authorised by **RESCO** Ltd.

1.4 PANEL GUARANTEE

Guarantee the compact laminate panels against mechanical failure of the resin panel for 10 years.

1.5 OTHER GUARANTEES

Hardware will carry the hardware manufacturer's normal guarantees.

Refer to the PRELIMINARIES AND GENERAL section for the required form of guarantee and details of when completed guarantees must be submitted.

2 PRODUCTS

2.1 PROPRIETARY LABORATORY BENCHTOP SYSTEM

There are 2 options the cost effective Multilab (mailab Bio) and the higher grade Laminex (Mailab Chem). Check the current literature to find out what bench is right for you!

3 EXECUTION

3.1 HANDLE AND STORE

Handle and store to avoid damage. Keep dry in transit. Store clear of and on a level floor and cover for protection.

3.2 PROTECT

Protect adjoining areas, surfaces and finish from damage.

3.3 INSTALL

Carry out the fixing, erection and fitting to finish rigid, plumb, true to line and face and square, to the requirements of RESCO NZ Ltd, according to the **RESCO** INSTALLATION INSTRUCTION MANUAL.

3.4 REPLACE

Replace any damaged or marked elements.

3.5 LEAVE

Leave work with parts fully and freely working and to the standard required by **RESCO** Ltd.

3.6 REMOVE

Remove debris, unused materials and elements from the site.